

Baby Cart

BC-100

Pure aluminum alloy column adopt the electrostatic spraying, the integral design is reasonable and artistic, convenient and durable.

Transparent and hi-intensity baby basin (PP Plastic). The baby basin adopts the healthy transparent resin material with environmental protection; it has good air permeability, healthy and environmental protection. The height of the baby crib is adjustable by the gas spring and the bed can be trendelenburged.

Security brakes.

Comes with smooth rotating casters for easy and silent transportation.

Inclination angle adjustable by gas spring, 0-12°
Height adjustable by gas spring, hi-low: 980mm-780mm.
Dimensions 870mm (L) 530mm (W) 940mm – 1090mm (H)



2 YEARS warranty

Neonatology

Success Through Quality/Since 1988

Advanced Instrumentations Inc. Success Through Quality, a Company You Can Trust

Advanced Instrumentations manufactures leading medical technology equipment in the areas of anesthesia, cardiology, operating room, gynecology and obstetrics, IV therapy, patient monitors, hospital furniture, neonatology and ultrasound. We deliver to the healthcare industry the highest-quality standards, reliability, and patient safety in all our products through effective, and rigorous testing procedures by our own department of Biomedical Engineering in the United States. All of our equipment comes with 2 years warranty and excellent post-sale support services.

Advanced Instrumentations Inc. Complies with the requirements of the ISO standards 9001: 2008 and 13485-2003 following the audit by one of the most prestigious global certification companies, as it is TÜV SÜD America. We comply with the requirements and are audited by the US Food and Drug Administration (FDA) an entity of the health and Human Services of the United States of America. These certifications are the result of dedication and commitment to excellence in our products and services.





6800 N.W. 77 Court, Miami, FL 33166 U.S.A.

Phone: 305-477-6331 Fax: 305-477-5351 2018 Advanced Instrumentations Inc., is a U.S.A registered company – All rights reserved.

All functionality, features, specifications and other product information provided in this document including, but not limited to, the benefits, design, pricing, components, performance, availability, and capabilities of the product are subject to change without notice or obligation. Advanced Instrumentations reserves the right to make changes to this document and the product described herein, at any time, without obligation on Advanced Instrumentations to provide notification of such change. Actual description and specification of the product in this document may be different. Images shown here are for representational purpose only, actual may vary.

Advanced and Advanced Instrumentations trademarks and logos shown are property of Advanced Instrumentations Inc.

