

# Baby Trolley

## BC-200

Instead of traditional mechanical buttons, adopt touch keyboard, which is more durable and sensitive.

The bed body can be multi-angle adjustment by switch, easy for mother and baby to communicate and feed.

Hand rail is easy to operate, multi-angle 360° adjustment. It's could be placed flat, or to angle 90°, can be used as a mosquito net support.

Tilt angle, 15°

Humanity design of Chassis, can be place diaper and other sundry.

4 pcs swivel castor, 2pcs with brake, 2pcs without.

Height adjustable by hydraulic systems, 855–1045mm.

Music playing: with a SD memory card, can storage music in MP3 format.

Weight scale: measure weight at anytime.

Atmosphere temperature measuring



**2 YEARS**  
warranty

Neonatology

2018 Advanced Instrumentations Inc, is a U.S.A registered company – All rights reserved.

All functionality, features, specifications and other product information provided in this document including, but not limited to, the benefits, design, pricing, components, performance, availability, and capabilities of the product are subject to change without notice or obligation. Advanced Instrumentations reserves the right to make changes to this document and the product described herein, at any time, without obligation on Advanced Instrumentations to provide notification of such change. Actual description and specification of the product in this document may be different. Images shown here are for representational purpose only, actual may vary.


Advanced and Advanced Instrumentations trademarks and logos shown are property of Advanced Instrumentations Inc.

## Success Through Quality/Since 1988

Advanced Instrumentations Inc.  
Success Through Quality,  
a Company You Can Trust

Advanced Instrumentations manufactures leading medical technology equipment in the areas of anesthesia, cardiology, operating room, gynecology and obstetrics, IV therapy, patient monitors, hospital furniture, neonatology and ultrasound. We deliver to the healthcare industry the highest-quality standards, reliability, and patient safety in all our products through effective, and rigorous testing procedures by our own department of Biomedical Engineering in the United States. All of our equipment comes with 2 years warranty and excellent post-sale support services.

Advanced Instrumentations Inc. Complies with the requirements of the ISO standards 9001: 2008 and 13485-2003 following the audit by one of the most prestigious global certification companies, as it is TÜV SÜD America. We comply with the requirements and are audited by the US Food and Drug Administration (FDA) an entity of the health and Human Services of the United States of America. These certifications are the result of dedication and commitment to excellence in our products and services.



6800 N.W. 77 Court,  
Miami, FL 33166  
U.S.A.

Phone: 305-477-6331  
Fax: 305-477-5351

For additional information visit us at:  
[www.advanced-inst.com](http://www.advanced-inst.com)