

Manual Operating Table

OT-30

OR



Operating Room

Operating Room


OT-30

Manual Operating Table

Universal Operating Tables for All Surgical Applications
The Advanced® OT-30 Operating Table is suitable for most surgical procedures including: chest, abdomen, nose, throat, gynecology and obstetrics, urology, orthopedics department. Advanced®OT-30 is used in Operating Rooms, Same Day Surgical Center and Labor & Delivery Center.

Features

- Excellent access for the surgical team to the patient.
- Versatile – use the same table for a wide variety of surgeries.
- Stainless steel construction with top quality materials for reliability and robustness.
- Easy to operate.
- Accessory rails, both sides.
- Adjustable height, hydraulic mechanism.
- Trendelenburg and reverse Trendelenburg adjustments .
- Lateral tilt adjustment (right).
- Flex inclination.
- Back section Up and Down.
- Leg section inclination Up Down.
- Leg plates horizontal opening.
- Adjustable head section.
- Central brake.
- One Table cushion can be removed for easy cleaning & disinfecting.
- One tabletop additional board, make easy to take X-Ray images.
- Meet ISO 13485 Quality Standard.
- Two years warranty.



2 YEARS
warranty

Technical Specifications

Physical Specifications	Dimensions Weight	2000mm Length X 500mm Width 293 Kg / 645 Lbs
Technical	Height without Pad Trendelenburg Reverse Trendelenburg Left Lateral Tilt Right Lateral Tilt Up Head Plate Down Head Plate Up Back Plate Down Back Plate Up Leg Plate Down Leg Plate Horizontal Leg Plate Flex Working Load	780mm -1020mm 27° 27° 21° 21° 60° 70° 75° 15° 15° 90° 90° 90° 225° 173Kg
Environmental	Temperature Relative Humidity Atmospheric Pressure	10~40°C 30%~75% 700~1060hpa
Storage Conditions	Temperature Relative Humidity Atmospheric Pressure	-40~70°C 10%~100% (including condensation) 500~1060 hpa
Configuration	Accessories Included	1 Anesthetic-screen 2 Shoulder support 2 Hand rest 1 Body support (long) 1 Body support (short) 2 Leg rest 1 Table cushion 2 Foot plate

2018 Advanced Instrumentations Inc, is a U.S.A registered company – All rights reserved.

All functionality, features, specifications and other product information provided in this document including, but not limited to, the benefits, design, pricing, components, performance, availability, and capabilities of the product are subject to change without notice or obligation. Advanced Instrumentations reserves the right to make changes to this document and the product described herein, at any time, without obligation on Advanced Instrumentations to provide notification of such change. Actual description and specification of the product in this document may be different. Images shown here are for representational purpose only, actual may vary.


Advanced and Advanced Instrumentations trademarks and logos shown are property of Advanced Instrumentations Inc.

Success Through Quality/Since 1988

Advanced Instrumentations Inc.
Success Through Quality,
a Company You Can Trust

Advanced Instrumentations manufactures leading medical technology equipment in the areas of anesthesia, cardiology, operating room, gynecology and obstetrics, IV therapy, patient monitors, hospital furniture, neonatology and ultrasound. We deliver to the healthcare industry the highest-quality standards, reliability, and patient safety in all our products through effective, and rigorous testing procedures by our own department of Biomedical Engineering in the United States. All of our equipment comes with 2 years warranty and excellent post-sale support services.

Advanced Instrumentations Inc. Complies with the requirements of the ISO standards 9001: 2008 and 13485-2003 following the audit by one of the most prestigious global certification companies, as it is TÜV SÜD America. We comply with the requirements and are audited by the US Food and Drug Administration (FDA) an entity of the health and Human Services of the United States of America. These certifications are the result of dedication and commitment to excellence in our products and services.



6800 N.W. 77 Court,
Miami, FL 33166
U.S.A.
Phone: 305-477-6331
Fax: 305-477-5351

For additional information visit us at:
www.advanced-inst.com