



Delivery Table OT-400

OR



Operating Room

Operating Room

OT-400

Delivery Table

The Advanced® OT-400 Delivery Table is used by OB/GYN professionals in Labor and Delivery Centers and other women health facilities .

Features

- Remote control/Easy to operate.
- Accessory rails, both sides.
- Power On/Off switch.
- Power On/Off status indicator light.
- Adjustable height.
- Trendelenburg / Reverse Trendelenburg adjustment.
- Electric central brake.
- Back section Up and Down.
- Removable leg section.
- Pad cover, easy cleaning and disinfection
- One Table cushion set, can be removed for easy cleaning & disinfecting.
- Power Supply: AC 100/240 Volts - 50/60 Hz.
- Meet ISO 13485 Quality Standard.
- Two years warranty.



2 YEARS
warranty

Technical Specifications

Physical Specifications	Dimensions Weight	1854mm X 570mm 330 Kg / 726Lbs
Technical	Height without Pad Trendelenburg Reverse Trendelenburg Up Back Plate Down Back Plate Working Load	762mm - 915mm 11° 25° 67° 13° 173 Kg
Power Supply	Voltage Frequency Power Capacity	110 V - 220 V 50Hz - 60Hz ≥1kw
Environmental	Temperature Relative Humidity Atmospheric Pressure	10°C~40°C ≤80% 760hPa~1060hPa
Storage Conditions	Temperature Relative Humidity Atmospheric Pressure	-30~55°C ≤95% 500~1060hPa
Configuration	Accessories Included	1 Remote Control 1 Table Cushion 2 Arm Support 2 Support Handles 2 Basic Stirrups 1 Leg Plate 1 Stainless Steel Drain Pan

2018 Advanced Instrumentations Inc, is a U.S.A registered company – All rights reserved.

All functionality, features, specifications and other product information provided in this document including, but not limited to, the benefits, design, pricing, components, performance, availability, and capabilities of the product are subject to change without notice or obligation. Advanced Instrumentations reserves the right to make changes to this document and the product described herein, at any time, without obligation on Advanced Instrumentations to provide notification of such change. Actual description and specification of the product in this document may be different. Images shown here are for representational purpose only, actual may vary.


Advanced and Advanced Instrumentations trademarks and logos shown are property of Advanced Instrumentations Inc.

Success Through Quality/Since 1988

Advanced Instrumentations Inc.
Success Through Quality,
a Company You Can Trust

Advanced Instrumentations manufactures leading medical technology equipment in the areas of anesthesia, cardiology, operating room, gynecology and obstetrics, IV therapy, patient monitors, hospital furniture, neonatology and ultrasound. We deliver to the healthcare industry the highest-quality standards, reliability, and patient safety in all our products through effective, and rigorous testing procedures by our own department of Biomedical Engineering in the United States. All of our equipment comes with 2 years warranty and excellent post-sale support services.

Advanced Instrumentations Inc. Complies with the requirements of the ISO standards 9001: 2008 and 13485-2003 following the audit by one of the most prestigious global certification companies, as it is TÜV SÜD America. We comply with the requirements and are audited by the US Food and Drug Administration (FDA) an entity of the health and Human Services of the United States of America. These certifications are the result of dedication and commitment to excellence in our products and services.



6800 N.W. 77 Court,
Miami, FL 33166
U.S.A.
Phone: 305-477-6331
Fax: 305-477-5351

For additional information visit us at:
www.advanced-inst.com