



Manual Stretcher

ST-100



FUR



Hospital Furniture

Hospital Furniture

ST-100

Manual Stetcher

Features:

Oxygen Cylinder holder

Four castors

Single-Direction retractable fifth wheel system Pedals in both sides

Stretcher pad

Foot-control brake system

One crank for bed up/down adjustment

Two IV Pole supports (head and foot)

One IV Pole

Pneumatic, lever-activated mechanism to lift/lower the back support

Powder coating paint

Meets ISO 13485 Quality Standard

Two years warranty



2 YEARS
warranty

Technical Specifications

Physical Specifications	Overall Dimensions	79.7" x 27.62" (2024.38mm x 701.54mm)
Specifications	Bedboard Dimensions	74" x 24" (1880mm x 610mm)
Specifications	Weight	99 Kg (217 Lb)
Adjustment	Hi - Low Back Rest Angle	22" ~ 33.7" (558.80mm ~ 856mm) 0° ~ 51°
	Side Rails	11" x 40.2" (279.40mm x 1016.50mm)
Wheels	Diameter	5.9" (150mm)
Weight	Capacity	250 Kg (550 Lb)

2018 Advanced Instrumentations Inc, is a U.S.A registered company – All rights reserved.

All functionality, features, specifications and other product information provided in this document including, but not limited to, the benefits, design, pricing, components, performance, availability, and capabilities of the product are subject to change without notice or obligation. Advanced Instrumentations reserves the right to make changes to this document and the product described herein, at any time, without obligation on Advanced Instrumentations to provide notification of such change. Actual description and specification of the product in this document may be different. Images shown here are for representational purpose only, actual may vary.


Advanced and Advanced Instrumentations trademarks and logos shown are property of Advanced Instrumentations Inc.

Success Through Quality/Since 1988

Advanced Instrumentations Inc.
Success Through Quality,
a Company You Can Trust

Advanced Instrumentations manufactures leading medical technology equipment in the areas of anesthesia, cardiology, operating room, gynecology and obstetrics, IV therapy, patient monitors, hospital furniture, neonatology and ultrasound. We deliver to the healthcare industry the highest-quality standards, reliability, and patient safety in all our products through effective, and rigorous testing procedures by our own department of Biomedical Engineering in the United States. All of our equipment comes with 2 years warranty and excellent post-sale support services.

Advanced Instrumentations Inc. complies with the requirements of the ISO standards 13585-2016 following the audit by one of the most prestigious global certification companies, as it is the SGS. We comply with the requirements and are audited by the US Food and Drug Administration (FDA) an entity of the Health and Human Services of the United States of America. These certifications are the result of dedication and commitment to excellence in our products and services.



6800 N.W. 77 Court,
Miami, FL 33166
U.S.A.

Phone: 305-477-6331
Fax: 305-477-5351

For additional information visit us at:
www.advanced-inst.com