

# Patient Stretcher

## ST-2000

FUR



Hospital Furniture

# Hospital Furniture

## ST-2000

### Patient Stetcher

The Advanced® ST-2000 hydraulic stretcher is designed for the efficient transport of patients in medical facilities, with height adjustable by pedals on both sides. The Advanced® ST-2000 hydraulic stretcher can carry patients weighing up to 550 pounds.

#### Features

- Steel frame
- X-Ray translucent plastic bed board
- X-Ray cassette tray ( optional )
- Stainless steel fold-down side rails
- Plastic chassis cover
- Polyurethane coated stretcher pad
- Oxygen Cylinder holder
- Four castors
- Single-Direction retractable fifth wheel system, pedals in both sides
- Central brake system with four brake pedals, in each corner
- Four IV Pole supports around the stretcher
- One IV Pole included
- Four rubber rolling bumper
- Five functions with control pedals in both sides
- Whole bed up/down
- Backrest up/down
- Leg part up/down
- Trendelenburg / Reverse Trendelenburg
- Pneumatic, lever-activated mechanism to lift/lower the back support
- One crank for adjusting knee position
- Pop-up push handles in the head side
- Large handle in the foot side
- Powder coating paint
- Meets ISO13485 Quality Standard
- Two years warranty



2 YEARS  
warranty

## Technical Specifications

Physical Specifications	Overall Dimensions	80.19" x 29.5" ( 2037mm x 749mm )
Specifications	Bedboard Dimensions	76.18" x 24.5" ( 1935mm x 622.3mm )
Adjustment	Hi-Low Trendelenburg Reverse Trendelenburg Back Rest Angle Knee Rest Angle	27" ~ 36.18" ( 686mm ~ 919mm ) ± 12° ± 12° 0° ~ 65° 0° ~ 45°
	Side Rails	17" x 43.5" ( 432mm x 1105mm )
Stretcher Pad	Dimensions Height	74.5" x 24.5" ( 1892.3mm x 622.3mm ) 4" ( 10cm )
Wheels	Diameter	7.8" ( 200mm )
Weight	Capacity	250 Kg ( 550 Lb )

2018 Advanced Instrumentations Inc, is a U.S.A registered company – All rights reserved.

All functionality, features, specifications and other product information provided in this document including, but not limited to, the benefits, design, pricing, components, performance, availability, and capabilities of the product are subject to change without notice or obligation. Advanced Instrumentations reserves the right to make changes to this document and the product described herein, at any time, without obligation on Advanced Instrumentations to provide notification of such change. Actual description and specification of the product in this document may be different. Images shown here are for representational purpose only, actual may vary.


Advanced and Advanced Instrumentations trademarks and logos shown are property of Advanced Instrumentations Inc.

## Success Through Quality/Since 1988

Advanced Instrumentations Inc.  
Success Through Quality,  
a Company You Can Trust

Advanced Instrumentations manufactures leading medical technology equipment in the areas of anesthesia, cardiology, operating room, gynecology and obstetrics, IV therapy, patient monitors, hospital furniture, neonatology and ultrasound. We deliver to the healthcare industry the highest-quality standards, reliability, and patient safety in all our products through effective, and rigorous testing procedures by our own department of Biomedical Engineering in the United States. All of our equipment comes with 2 years warranty and excellent post-sale support services.

Advanced Instrumentations Inc. complies with the requirements of the ISO standards 13585-2016 following the audit by one of the most prestigious global certification companies, as it is the SGS. We comply with the requirements and are audited by the US Food and Drug Administration (FDA) an entity of the Health and Human Services of the United States of America. These certifications are the result of dedication and commitment to excellence in our products and services.



6800 N.W. 77 Court,  
Miami, FL 33166  
U.S.A.

Phone: 305-477-6331  
Fax: 305-477-5351

For additional information visit us at:  
[www.advanced-inst.com](http://www.advanced-inst.com)