# Saturn 8000 17x17

Flat Panel Detector for Digital Radiography





Saturn 8000 17x17 is a NMI's flat panel detector with a large field coverage area of 17"x17", which is designed for general radiographic application using its unique image processing system and proprietary flat panel detector. The active 9.4 Mega-pixel sensor ensures superior image quality to meet market demand for precise diagnosis performance, which is made possible by the NMI's digital radiographic flat panel detector technology. With automatic exposure detection, Saturn 8000 17x17 can be used without connection to the X-ray generator. It works just like CR or film as it is triggered by the radiation, instead of a generator interface connection.



## Saturn 8000 17x17

Flat Panel Detector for Digital Radiography



#### **Features**

- \* Wide active area of 17" x 17"
- \* High spatial resolution with 140um pixel array
- \* Stable and reliable automatic exposure detection
- \* Viewer software running on Windows™ OS (DxWorks)
- \* Communication interface through Gigabit Ethernet (1000 BASE-T)
- \* Simple and easy integration with all kinds of digital radiography systems



Configuration Detector Power Supply Unit

### **Acquired Images**





### Drawing

470





## **Technical Specification**

Application General radiography

Technology Flat panel detector: a-Si TFT with PIN diode

Scintillator Csl:Tl / Gd,O,S:Tb Pixel Pitch 140um x 140um Pixels 3,072 x 3,072 pixels Image Size 17 x 17 inches (43 x 43cm)

A/D Conversion 14 bit Grayscale 16,384 steps X-ray Voltage Range 40 ~ 150kVp

X-ray Generator Interface Line trigger: DR Trigger Mode

Auto trigger: AED (Automatic Exposure Detection) Mode

Data Interface Gigabit Ethernet (1000BASE-T) Dimensions 470 (W) x 470 (L) x 35 (T) mm

Weight Approx. 11kg

Operating Environment 15 ~ 35 °C, 30 ~ 80% RH (non-condensing)

Power DC24V, 0.8A

\* Specifications are subject to change without prior notification.







#### Office

3F., NO.32, Sec. 2, Zhonyang S. Rd., Beitou Dist., Taipei City 11270, Taiwan Tel +886-2-2892-2275 Fax +886-2-2892-2172 E-mail Sales@newmedical.com.tw

#### Factory

41-3, Burim-ro 170beon-gil, Dongan-gu, Anyang-si, Gyeonggi-do, 14055 Republic of KOREA

www.nmi-dr.com

