



GP5 – 4 WAY FLOAT TOP TABLE SYSTEM WITH ROTATE TUBE STAND

Choice of 20, 32, 40 or 50 kW High Frequency Generator

HF Generators


SHF210AP – 320 mA, 125 kVp (250 mA @ 80 kVp), 16 mA Stations
 SHF310AP – 400 mA, 125 kVp (400 mA @ 80 kVp), 17 mA Stations
 SHF410AP – 500 mA, 125 kVp (500 mA @ 80 kVp), 18 mA Stations
 SHF510AP – 640 mA, 125 kVp (640 mA @ 78 kVp), 19 mA Stations

- 0.001 to 10 Second Time Range
- 0.1 to 500 mAs
- kVp Selection, 40 to 125 kVp in 1 kV Increments
- Digital Displays
- Automatic Line Compensation
- APR, Anatomical Programming
- Table Mount
- Exposure Hand Switch
- UL/CE/FDA/Listed
- For Operation on 1 Phase 208/240VAC – 3 Phase Optional



POWER CABINET SPECIFICATIONS

Line Powered



MULTI-FWFTT-G Four Way Float Top Table, Fixed Height



- 86.5" x 32.4" (2200 mm x 825 mm) Float Top Table
- Table Top Supports Patients up to 440 lbs (200 kgs)
- Longitudinal Table Top Travel 35" (90 cm), $\pm 17.5"$
- Transverse Table Top Travel 9" (23 cm), $\pm 4.5"$
- Table Top to Film Distance 30.3" (76.9 cm)
- Electromagnetic Brakes
- Heavy Duty Grid Cabinet (Bucky Optional) with 103 Line, 10:1 Ratio Grid and Cassette Tray

13-CSM-84 Rotate Tube Stand

- Floor-to-Wall Mounted
- $\pm 5"$ Transverse Travel, 180 Column Rotation
- Vertical Travel 14.5" to 76" (36.8 to 193 cm)
- Fits in Ceiling Heights of 97" to 102" (246.4 cm to 261.6 cm)
- Platform Tube Mount, Angulation Dial and Operator Handles
- Magnetic Locks Control Longitudinal, Transverse, Vertical and Tube Rotation
- 80" (203.2 cm) Longitudinal Travel

13-BF-7E Wall Stand with Grid Cabinet

- 84" (213.4 cm) Floor-to-Wall Mounted Column
- Magnetic Locks
- 17 x 17" (43 x 43 cm) Heavy Duty Grid Cabinet with 103 Line, 10:1 Ratio Grid and Cassette Tray
- Bucky, Various Grids and Tilting Mechanism Optional



E-7239FX X-Ray Tube

- 140,000 Heat Unit
- 1.0 x 2.0 mm Focal Spots

CM-150-24-C125 Collimator

- Manual Collimator
- Laser Light
- Swivel Mount

10540-25 High Voltage Cables

- 25' (762 cm) Long

TXR • 5481 Skyland Blvd., E • Cottondale, AL 35453

Phone: 205-556-3803 • Fax: 205-556-3824 • 800-897-9729 • www.txr.com • sales@txr.com

Complies with FDA and Product Safety Requirements. • Specifications Subject to Change Without Notice

0415-GP5